











# Complete Laboratory Catalogue



# Turning Technology into Solutions



# Precise Results Unrivalled in its Performance Class

No matter you are carrying out a simple routine task or a complicated weighing procedure, New Citizen Balances will impress you with their excellent performance, delivering accurate results within the shortest possible time, thanks to MonoTech and Perfect Self Automatic Calibration (PSAC).

The proven MonoTech weighing technology guarantees long operation life and extreme ruggedness, while PSAC automatically calibrate the balances using internal weights whenever changes in temperature affect the accuracy of your weighing results. For added security, you can configure PSAC to initiate self-calibration at specific time and every power on.

## Built-in application programme other than basic weighing ease your routine lab work like,

**%** 

Percent Weighing



Piece Counting



Checkweighing



**Unit Conversion** 



**Animal Weighing** 



Statistic



Totalization



Formulation



Density Determination





#### **EXCELLENT READABILITY**

High Contrast Backlit LCD display is exceptionally easy to read under any room lighting conditions.

20-Jul-10	10:32AM
Citizen Model	Scales CX 220
Ser.no.	9223102
Ver.no.	r0.1.5.3
ID	1234567
Calibration:	Internal
Start:	Manua1
Temperature	29.449'C
Diff.	+ 0.0009g
Internal Cal	
Diff.	0.0000g
20-Jul-10 Name:	10:32AM

## STATISTICS:

This smart feature enables the user to obtain the statistics of the data stored in the scale including details such as number of readings, minimum value, maximum value, standard deviation, difference and Coefficient of variance



#### THE BRILLIANT MONO-TECH

Unique & unequalled high end MONO-TECH technology (Weigh Sensor) to match the highest international standards.



# METAL HOUSING

The metal housing protects the weighing sensor from environmental influences which enhances the performance even with intensive use in harsh environments



## ELECTROMAGNETIC COMPATIBILITY

Metal housing eliminates the effects of thunderstorms or electromagnetic inter-ference from other equipments.



# **EASY DATA TRANSFER**

Only a single cable is required to transfer balance data to MSExcel or other Windows application.

No troublesome interface attachment or software installation is necessary. This unique feature makes Citizen balance extremely convenient for data analysis.

Optional USB Interface also available on request at the time



r	
28-Jul-10	03:19PM
Citizen Model Ser.no. Ver.no. ID	CX 220 9223102 r0.1.5.3 1234567
1 + 2 + 3 + 4 +	1.5750 g 2.3500 g 7.8950 g 4.1750 g
n max min avg std var	4 7.8950 g 1.5750 g 3.9987 2.8169 7.9349
28-Ju1-10 Name:	03:23PM

# ISO/GLP Compliant Printouts:

Built-in real time clock Makes sure that the data are the true reflection of the test results & the data is traceable document in compliance to SOP's and GLP regulation for use in quality management systems.



#### **Carat Balance**

Has been designed to meet high requirement for measuring precious gems, with high readability of 1 mili carat.



## **Density Balance**

Has been designed to measure density / specific gravity of solid material (Purity of Gold & Platinum in Karat).

# Weigh Below



Citizen's built-in weigh-below hook allows you to suspend an item below the balance to calculate density / specific gravity.



Large draft shield chamber, gives more space to the user, samples can be easily placed on the pan & also it gets accustomed to the interior temperature.

Optimal lighting conditions inside the draft shield, thanks to its nearly frame less all-glass design.





	MODEL	CAPACITY	READABILITY	REPEATABILITY (SD)	LINEARITY	PAN SIZE (mm)	RESPONSE TIME	CAL
Analytical Balance	CX 301	301 g	0.1 mg	0.2 mg	0.3 mg	90 <b>Ø</b>	3 sec.	Internal
2.7	CX 220	220 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	Internal
	CX 120	120 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	Internal
	CX 65	65 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	Internal
	CY 204	220 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	External
	CY 104	120 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	External
DIMEN: 345 x 212 x 341	CY 64	60 g	0.1 mg	0.1 mg	0.2 mg	90 <b>Ø</b>	3 sec.	External
Precision Balance	CY 120	120 g	0.001g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
	CY 220	220 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
	CY 320	320 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
Comment of the second of the s	CY 420	420 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
	CY 510	510 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
	CY 720	720 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	External
	CY 1003	1000 g	0.001 g	0.002 g	0.002 g	128 x 128	2-3 sec.	External
	CY 220C	220 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	Internal
	CY 320C	320 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	Internal
DIMEN: 342.5 x 212 x 271	CY 510C	510 g	0.001 g	0.001 g	0.002 g	128 x 128	2-3 sec.	Internal
	CY 1003C	1000 g	0.001 g	0.002 g	0.002 g	128 x 128	2-3 sec.	Internal\
Top Loading Precision Balance	CG 1202	1200 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	External
COLISON COLINE	CG 2202	2200 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	External
	CG 3102	3100 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	External
	CG 4102	4100 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	External
	CG 6102	6100 g	0.01 g	0.02 g	0.03 g	190 x 205	2-3 sec.	External
	CG 2202C	<b>220</b> 0 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	Internal
	CG 4102C	4100 g	0.01 g	0.01 g	0.02 g	190 x 205	2-3 sec.	Internal
DIMEN: 342.5 x 212 x 89.5	CG 6102C	6100 g	0.01 g	0.02 g	0.03 g	190 x 205	2-3 sec.	Internal
<b>High Capacity Precision Balance</b>	SSH 95	15 Kg	0.1g	0.1 g	0.2 g	400 x 300	2-3 sec.	External
	SSH 94	20 Kg	0.1g	0.1 g	0.2 g	400 x 300	2-3 sec.	External
	SSH 94H	25 Kg	0.1g	0.1 g	0.2 g	400 x 300	2-3 sec.	External
	SSH 93	35 Kg	0.1g	0.1 g	0.2 g	400 x 300	2-3 sec.	External

1g

1g

1 g

1 g

2 g

2 g

400 x 300

400 x 300

2-3 sec.

2-3 sec.

External

External





DIMEN: 517 x 302 x 130

#### WORLD WIDE MARKETING

# CITIZEN SCALE (I) PVT. LTD.

SSH 92

SSH 91

50 Kg

60 Kg

Citizen House, Unit No. E-2, Plot No. 11, WICEL,

Opp. Seepz Gate No. 1, Andheri (East), Mumbai - 400 093. (INDIA)

Tel.: +91 - 22 - 4243 7700 | Fax: +91 - 22 - 4243 7800

E-mail: sales@citizescales.com | Web: www.citizenscales.com

